Troubleshooting > Atellica CH > Troubleshooting Instructions > Service Software

SIEMENS																																																																																																																																																																																																																																															
Home Health Instrument Autocheck Autocheck Subsystem Subsystem Auto Align Detectors Prime Sequences Devices Charts Info CanNetwork Thermals WebMon DB-Tools																																																																																																																																																																																																																																															
Start Autocheck											Rectangular Strip																																																																																																																																																																																																																																				
Id: GCCEP2 Type: CC Model: CC1 Hardware lot 5 125.3.4425000 - (kachro20 - GGW7EKACHERIES) 2022-06-02 15:16:05																																																																																																																																																																																																																																															
<table border="1"> <thead> <tr> <th>Command test and Read device status</th><th>Number</th><th>Read Pass</th><th>Read Fail</th><th>Hardware Emulation</th><th>Complementary device tests</th><th>Number</th><th>Pass</th><th>Fail</th><th>Skipped</th><th></th><th></th></tr> </thead> <tbody> <tr> <td>Digital Inputs</td><td>32</td><td>32</td><td>0</td><td>0</td><td>Alignments</td><td>105</td><td>105</td><td>0</td><td>0</td><td></td><td></td></tr> <tr> <td>Digital Outputs</td><td>82</td><td>82</td><td>0</td><td>0</td><td>Water Sensors</td><td>12</td><td>11</td><td>1</td><td>0</td><td></td><td></td></tr> <tr> <td>Analog Inputs</td><td>70</td><td>69</td><td>1</td><td>0</td><td>Probe Sensors</td><td>10</td><td>10</td><td>0</td><td>0</td><td></td><td></td></tr> <tr> <td>Motors</td><td>32</td><td>32</td><td>0</td><td>0</td><td>Bath Level Sensor</td><td>1</td><td>1</td><td>0</td><td>0</td><td></td><td></td></tr> <tr> <td>Boards</td><td>23</td><td>23</td><td>0</td><td>0</td><td>Lowest Working Current +</td><td>11</td><td>11</td><td>0</td><td>0</td><td></td><td></td></tr> <tr> <td>Connectivity and drive test</td><td>Off Test Pass</td><td>Test Pass</td><td>On Test Pass</td><td>Fail</td><td>Motor Backlash +</td><td>24</td><td>24</td><td>0</td><td>0</td><td></td><td></td></tr> <tr> <td>Digital Outputs (DC)</td><td>57</td><td>0</td><td>57</td><td>0</td><td>Transmission Check</td><td>186</td><td>186</td><td>0</td><td>0</td><td></td><td></td></tr> <tr> <td>Thermals</td><td>8</td><td>0</td><td>8</td><td>0</td><td>Probe Rinse Pressure +</td><td>4</td><td>4</td><td>0</td><td>0</td><td></td><td></td></tr> <tr> <td>Analog Inputs (Nominal)</td><td></td><td></td><td>32</td><td>0</td><td>Probe Level Sense +</td><td>14</td><td>12</td><td>2</td><td>0</td><td></td><td></td></tr> <tr> <td>System Fluid Valves</td><td></td><td></td><td>4</td><td>0</td><td>Probe Alignment Check</td><td>4</td><td>4</td><td>0</td><td>0</td><td></td><td></td></tr> <tr> <td>Utilities</td><td>Test</td><td>Value</td><td></td><td></td><td>Probe leak test +</td><td>4</td><td>4</td><td>0</td><td>0</td><td></td><td></td></tr> <tr> <td>Vacuum Pump +</td><td>Pass</td><td>6.65 inHg</td><td>3087</td><td></td><td>Dynamic Fluid Pressure +</td><td>5</td><td>5</td><td>0</td><td>0</td><td></td><td></td></tr> <tr> <td>Fixed Orifice Rise +</td><td>Pass</td><td></td><td></td><td></td><td>Analog pairs</td><td>16</td><td>13</td><td>1</td><td>2</td><td></td><td></td></tr> <tr> <td>UPS +</td><td>Pass</td><td></td><td></td><td></td><td>Saline Pump test +</td><td>1</td><td>1</td><td>0</td><td>0</td><td></td><td></td></tr> <tr> <td>Database Size +</td><td>Pass</td><td></td><td></td><td></td><td>Bleach Pump test +</td><td>1</td><td>1</td><td>0</td><td>0</td><td></td><td></td></tr> <tr> <td>Detectors</td><td>Tests</td><td>Pass</td><td>Fail</td><td>Skipped</td><td>Dilution Washer Testing +</td><td>11</td><td>11</td><td>0</td><td>0</td><td></td><td></td></tr> <tr> <td>Photometer +</td><td>44</td><td>43</td><td>1</td><td>0</td><td>Reaction Washer Testing +</td><td>19</td><td>19</td><td>0</td><td>0</td><td></td><td></td></tr> <tr> <td>IMT +</td><td>4</td><td>3</td><td>1</td><td>0</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr> </tbody> </table>												Command test and Read device status	Number	Read Pass	Read Fail	Hardware Emulation	Complementary device tests	Number	Pass	Fail	Skipped			Digital Inputs	32	32	0	0	Alignments	105	105	0	0			Digital Outputs	82	82	0	0	Water Sensors	12	11	1	0			Analog Inputs	70	69	1	0	Probe Sensors	10	10	0	0			Motors	32	32	0	0	Bath Level Sensor	1	1	0	0			Boards	23	23	0	0	Lowest Working Current +	11	11	0	0			Connectivity and drive test	Off Test Pass	Test Pass	On Test Pass	Fail	Motor Backlash +	24	24	0	0			Digital Outputs (DC)	57	0	57	0	Transmission Check	186	186	0	0			Thermals	8	0	8	0	Probe Rinse Pressure +	4	4	0	0			Analog Inputs (Nominal)			32	0	Probe Level Sense +	14	12	2	0			System Fluid Valves			4	0	Probe Alignment Check	4	4	0	0			Utilities	Test	Value			Probe leak test +	4	4	0	0			Vacuum Pump +	Pass	6.65 inHg	3087		Dynamic Fluid Pressure +	5	5	0	0			Fixed Orifice Rise +	Pass				Analog pairs	16	13	1	2			UPS +	Pass				Saline Pump test +	1	1	0	0			Database Size +	Pass				Bleach Pump test +	1	1	0	0			Detectors	Tests	Pass	Fail	Skipped	Dilution Washer Testing +	11	11	0	0			Photometer +	44	43	1	0	Reaction Washer Testing +	19	19	0	0			IMT +	4	3	1	0							
Command test and Read device status	Number	Read Pass	Read Fail	Hardware Emulation	Complementary device tests	Number	Pass	Fail	Skipped																																																																																																																																																																																																																																						
Digital Inputs	32	32	0	0	Alignments	105	105	0	0																																																																																																																																																																																																																																						
Digital Outputs	82	82	0	0	Water Sensors	12	11	1	0																																																																																																																																																																																																																																						
Analog Inputs	70	69	1	0	Probe Sensors	10	10	0	0																																																																																																																																																																																																																																						
Motors	32	32	0	0	Bath Level Sensor	1	1	0	0																																																																																																																																																																																																																																						
Boards	23	23	0	0	Lowest Working Current +	11	11	0	0																																																																																																																																																																																																																																						
Connectivity and drive test	Off Test Pass	Test Pass	On Test Pass	Fail	Motor Backlash +	24	24	0	0																																																																																																																																																																																																																																						
Digital Outputs (DC)	57	0	57	0	Transmission Check	186	186	0	0																																																																																																																																																																																																																																						
Thermals	8	0	8	0	Probe Rinse Pressure +	4	4	0	0																																																																																																																																																																																																																																						
Analog Inputs (Nominal)			32	0	Probe Level Sense +	14	12	2	0																																																																																																																																																																																																																																						
System Fluid Valves			4	0	Probe Alignment Check	4	4	0	0																																																																																																																																																																																																																																						
Utilities	Test	Value			Probe leak test +	4	4	0	0																																																																																																																																																																																																																																						
Vacuum Pump +	Pass	6.65 inHg	3087		Dynamic Fluid Pressure +	5	5	0	0																																																																																																																																																																																																																																						
Fixed Orifice Rise +	Pass				Analog pairs	16	13	1	2																																																																																																																																																																																																																																						
UPS +	Pass				Saline Pump test +	1	1	0	0																																																																																																																																																																																																																																						
Database Size +	Pass				Bleach Pump test +	1	1	0	0																																																																																																																																																																																																																																						
Detectors	Tests	Pass	Fail	Skipped	Dilution Washer Testing +	11	11	0	0																																																																																																																																																																																																																																						
Photometer +	44	43	1	0	Reaction Washer Testing +	19	19	0	0																																																																																																																																																																																																																																						
IMT +	4	3	1	0																																																																																																																																																																																																																																											
Complete 10:14 Total Tests: 292 Total Test Time: 58:49																																																																																																																																																																																																																																															

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
Figure 4: High Pressure Transducer Failed


Figure 5: Water Bottle Empty

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Tab: Digital Outputs

Subsystem	Device	Value
DilutionArm	probe bleach pump	650
DilutionArm	probe drain valve	186
DilutionArm	probe imt drain valve	186
DilutionArm	probe inlet valve	187
DilutionArm	probe pump valve	189
DilutionArm	saline pump	58
DilutionArm	saline valve	192
DilutionMixer	probe drain valve	187
DilutionWasher	washer 1 solution pump	477
DilutionWasher	washer 1 waste valve	185
DilutionWasher	washer 1 water valve	187
DilutionWasher	washer 2 waste valve	185
DilutionWasher	washer 2 water valve	186
DilutionWasher	washer 3 waste valve	185
DilutionWasher	washer 3 water valve	186
IMT	flush pump std A	643
IMT	imt diluent valve	193
IMT	salt bridge valve	191
LiquidWaste	waste pump	457
LysingAgentThermal	ted fan	139
ReactionRingBathLevel	drain valve	261
ReactionRingBathLevel	water valve	190
ReactionWasher	washer 1 waste valve	181
ReactionWasher	washer 1 water valve	190
ReactionWasher	washer 2 solution pump	474
ReactionWasher	washer 2 waste valve	184
ReactionWasher	washer 2 water valve	191
ReactionWasher	washer 3 solution pump	474
ReactionWasher	washer 3 waste valve	181
ReactionWasher	washer 3 water valve	190
ReactionWasher	washer 4 waste valve	181
ReactionWasher	washer 4 water valve	189
ReactionWasher	washer 5 waste valve	181
ReactionWasher	washer 5 water valve	190
ReagentArm1	probe drain valve	186
ReagentArm1	probe inlet valve	187
ReagentArm1	probe pump valve	191
ReagentArm2	probe drain valve	186

Subsystem	Device	Value
ReagentArm2	probe inlet valve	187
ReagentArm2	probe pump valve	191
ReagentLoaderArm	gripper washer drain valve	186
ReagentLoaderArm	gripper washer water valve	186
ReagentMixer	probe drain valve	186
ReagentServer1	ted 1 fan	387
ReagentServer1	ted 2 fan	399
ReagentServer2	ted 1 fan	375
ReagentServer2	ted 2 fan	386
SampleArm1	probe drain valve	186
SampleArm1	probe inlet valve	187
SampleArm1	probe pump valve	190
SampleMixer	probe drain valve	187
SystemFluidCuvetteConditioner	bottle valve	186
SystemFluidCuvetteWash	bottle valve	187
SystemFluidProbeCleaner	bottle valve	187
SystemFluidSalineDiluent	bottle valve	187
WaterSupply	degasser vacuum pump	468
WaterSupply	water supply inlet valve	290

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Subsystem	Device	Value
MixerThermal	thermal fluid line verify thermistor	33
ProbeThermal	thermal fluid line control thermistor	32
ProbeThermal	thermal fluid line verify thermistor	32
ReactionRingBathThermal	reaction ring control thermistor	37
ReactionRingBathThermal	reaction ring verify thermistor	37
ReagentArm1	pressure transducer (HSen)	8981
ReagentArm1	pressure transducer (LSen)	1707
ReagentArm2	pressure transducer (HSen)	8962
ReagentArm2	pressure transducer (LSen)	1703
ReagentServer1Thermal	reagent server ted 1 thermistor	8
ReagentServer1Thermal	reagent server ted 2 thermistor	9
ReagentServer2Thermal	reagent server ted 1 thermistor	9
ReagentServer2Thermal	reagent server ted 2 thermistor	9
SampleArm1	pressure transducer (HSen)	9019
SampleArm1	pressure transducer (LSen)	1713
Vacuum	manifold pressure transducer	5345
WaterSupply	high pressure transducer	10
WaterSupply	inlet monitor thermistor	29
WaterSupply	inlet verify thermistor	29
WaterSupply	low pressure transducer	10
WaterSupply	recirculation monitor thermistor	32
WaterSupply	recirculation verify thermistor	32
WaterSupply	tank level sensor	82

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Autocheck Subsystem									
Start Autocheck									
Id: GGCCEP2 Type: CC Model: CC1 Hardware lot: 5 1.22.0.3637000 - (kachro20 - OGW7EKACHEIRIES) 2020-03-23 10:21:21									
<input type="checkbox"/> Loop while checked									
Command test and Read device status									
Number		Read Pass	Read Fail	Hardware Emulation	Complementary device tests		Number	Pass	Fail
Digital Inputs	32	32	0	0	Alignments	89	89	0	0
Digital Outputs	82	82	0	0	Water Sensors	12	12	0	0
Analog Inputs	70	70	0	0	Probe Sensors	10	10	0	0
Motors	32	32	0	0	Bath Level Sensor	1	1	0	0
Boards	23	23	0	0	Lowest Working Current +		11	10	1
Connectivity and drive test									
Off Test Pass		On Test Pass	On Test Fail	Skipped	Motor Backlash +		24	24	0
Digital Outputs (DC)	57	0	57	0	Transmission Check		186	186	0
Thermals	9	0	9	0	Probe Rinse Pressure +		4	4	0
Analog Inputs (Nominal)			33	0	Probe Level Sense +		5	5	0
Utilities	Test		Value		Probe Alignment Check		4	4	0
Vacuum Pump +	Pass		6.75 inHg 3132		Probe leak test +		4	4	0
Pressure Relief Valve +	Pass				Dynamic Fluid Pressure +		5	5	0
UPS +	Pass				Analog pairs		16	15	1
Database Check +	Pass				Saline Pump test +		1	1	0
Detectors									
Photometer +	Tests	Pass	Fail	Skipped	Bleach Pump test +		1	1	0
IMT +	43	43	0	0	Dilution Washer Testing +		11	11	0
	4	3	1	0	Reaction Washer Testing +		19	19	0

Figure 6: Autocheck - Vacuum Pump

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
Figure 7: Vacuum Pump Test Detail Screen

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
Figure 8: Autocheck Tab - Results - Pressure Relief Valve

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
Figure 9: Pressure Relief Valve Graphing - Good vs. Bad

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- **Battery Remaining Capacity (%)**
- **Battery Run Time (Minutes)**
- **Battery Voltage (V)**

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View Autocheck UPS Results		
Manufacturer	Emerson Network Power	
Product	Liebert GXT4	
Serial Number	1429000256AGBE3	
Firmware	130:110:0	
Parameter	Input	Output
Voltage (V)	225	231
Current (A)	0.600	0.400
Frequency (Hz)	60	60
Config Voltage (V)		230
Load (%)		14
Test	Before	After
Battery Remaining Capacity (%)	100	99
Battery Run Time (minutes)	51	46
Battery Voltage (V)	54	50
Battery Test		Pass

Figure 10: Autocheck Tab - Results - UPS

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Concentration = pump volume / water volume

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water volume = flow rate \times valve on time

View Autocheck Cuvette Washer Test Results									
	Min	Max	Probe 1	Probe 2	Probe 3	Probe 4	Probe 5	Probe 6	Dryer
Full (µl)			460	460	460	460	460	460	230
Skim (µl)	300	425	372	366	337	365	369
Wash (µl)	140	415	304	307	347	367	387
Empty (µl)		100	43	14	29	14	16	15	3
Pump (µl)	26	60	...	37.8
Pump (µl)	7	30	19.7
Concentration (1: X)			...	9.4	19.6
Vacuum Ratio	0.320	0.760	0.520
Valve On time (ms)			248	315	360	359	359
Ideal On time (ms)			324	361	329	295	292

Figure 36: Autocheck Tab - Results - Reaction Washer Testing

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Figure 37: Reaction Washer Results - High Altitude (>2000 m)

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Overview, **Autochecks** и **Health Reports**.

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The screenshot shows a table with two main sections. The left section is a summary table with columns: Date / Time, Launch, Tests, and Result. The right section is a detailed log table with columns: Result, Test Name, Version, Device, Sensor, Data Name, and Value. The log table contains numerous rows of data, mostly 'FAIL' entries, with some 'PASS' and 'ERROR' entries interspersed.

Date / Time	Launch	Tests	Result			
1/7/2017 4:31 AM	PROC_AUTO_S	1	ERROR			
Result	Test Name	Version	Device	Sensor	Data Name	Value
FAIL	Level Sense	1	DilutionArm:Drain		Width	92
FAIL	Level Sense	1	DilutionArm:IMTDrain		Width	92
FAIL	Level Sense	1	ReagentArm1:Drain		Width	92
FAIL	Level Sense	1	ReagentArm2:Drain		Width	92
FAIL	Level Sense	1	SampleArm1:Drain		Width	92
FAIL	Level Sense	1	DilutionArm:Drain		Delay	1001.2
FAIL	Level Sense	1	DilutionArm:IMTDrain		Delay	1005
FAIL	Level Sense	1	ReagentArm1:Drain		Delay	1001.2
FAIL	Level Sense	1	ReagentArm2:Drain		Delay	1000
FAIL	Level Sense	1	SampleArm1:Drain		Delay	1001.2
FAIL	Probe Rinse Press...	1	ReagentArm1		Pressure drop	0.016
FAIL	Probe Rinse Press...	1	SampleArm1		Pressure drop	0.015
FAIL	Probe Rinse Press...	1	ReagentArm2		Pressure drop	0.014
FAIL	Probe Rinse Press...	1	DilutionArm		Pressure drop	-0.017
FAIL	Probe Rinse Press...	1	ReagentArm1		Pressure	0.01
FAIL	Probe Rinse Press...	1	SampleArm1		Pressure	0.02
FAIL	Probe Rinse Press...	1	ReagentArm2		Pressure	0.04
FAIL	Probe Rinse Press...	1	DilutionArm		Pressure	0
FAIL	Pump	1	DilutionArm:probe b...		Volume	0
FAIL	Pump	1	DilutionArm:saline p...		Ratio	0
PASS	Autocheck	1	Elapsed Time			331.26
PASS	Computer	1	VMware, Inc.		VMware Virtu...	3
PASS	Cuvette Status	1	ReactionRing		Number Of T...	221
PASS	Cuvette Status	1	ReactionRing		Number Of Q...	0
PASS	Cuvette Status	1	ReactionRing		Number Of N...	0
PASS	Database Check	1	Module		Errors	0

Figure 40: Autocheck Results

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The screenshot shows a table with two main sections. The left section is a summary table with columns: Created, Start, End, Launch, View, and Export. The right section is a form titled 'Create New Health Report' with fields for Start Time and End Time, and a 'Create New' button.

Created	Start	End	Launch	View	Export
1/7/2017 4:06 AM	12/31/2016 4:06 AM	1/7/2017 4:06 AM	Maintenance	View	Export

Create New Health Report

Start Time: End Time:

Create New

Figure 41: Health Tab - Health Reports Sub-tab

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
Figure 42: Health Tab - Health Reports Sub-tab - Report Complete

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
Figure 43: Health Tab- Health Reports Sub-tab - View / Export Report

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
Figure 44: Health Tab - Health Reports Sub-tab - Export Report Floppy Message

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
Figure 45: Health Tab - Health Reports Sub-tab - Export Location

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